

## Reax 2 con adaptador universal

### Technical datasheet



|                                               |                                             |
|-----------------------------------------------|---------------------------------------------|
| <b>Dimensions (W × H × D)</b>                 | 501 × 235 × 180 mm                          |
| <b>Weight</b>                                 | approx. 5 kg                                |
| <b>Maximum permissible load</b>               | 1 kg                                        |
| <b>Drive</b>                                  | Shaded pole motor                           |
| <b>Speed range</b>                            | 20 – 200 rpm; Indicator: analogue           |
| <b>Protection class (EN 60529)</b>            | IP21                                        |
| <b>Acoustic pressure</b>                      | < 45 (dB(A)) (in accordance with IEC 61010) |
| <b>Rated voltage</b>                          | 230-240 V (50 Hz) or 115 V (60 Hz)          |
| <b>Protection class</b>                       | I                                           |
| <b>Overvoltage category</b>                   | n.a.                                        |
| <b>Degree of pollution</b>                    | n.a.                                        |
| <b>Power input</b>                            | 33 W max.                                   |
| <b>Operating temperature</b>                  | n.a.                                        |
| <b>Installation altitude</b>                  | n.a.                                        |
| <b>Overheat protection</b>                    | self reset function                         |
| <b>Speed control</b>                          | electronic                                  |
| <b>Usable area</b>                            | 230 × 160 mm                                |
| <b>Movement</b>                               | Overhead rotation                           |
| <b>Orbit</b>                                  | n.a.                                        |
| <b>Fuse</b>                                   | n.a.                                        |
| <b>Permissible mains voltage fluctuations</b> | n.a.                                        |

**Bottle dimensions**

n.a.